

Menicon PremiO

Silicone hydrogel contact lenses

2 week
replacement

Material

MeniSilk™ & Nanogloss™ Technology

- Asmofilcon A.
- Silicone Hydrogel (Group V ISO18369-1).
- Water content: 40%.
- **Dk: 129×10^{-11}** (cm²/sec) (mLO₂/(mL · mmHg)) (ISO).
- Spectral and luminous transmittance: 98%.
- Refractive index: 1.423.
- Contact angle: 27° (Captive bubble method).
- Blue handling tint.
- Nanogloss™ surface treatment.

Parameters range

	premiO		premiO ^{toric}		
Ø _r	14.00 mm		14.00 mm		
r ₀	8.60 mm	8.30 mm	8.60 mm		
t _c	0.08 mm (-3.00D)		0.08 mm (-3.00D)		
Dk/t _c	161 x 10⁻⁹ (cm/sec) • (mLO ₂ /(mL x mmHg))		161 x 10⁻⁹ (cm/sec) • (mLO ₂ /(mL x mmHg))		
Sph P	+6.00 ~ - 6.25 D (in 0.25) -6.50 ~ - 13.00D (in 0.50)		plan ~ -6.00D (in 0.25) -6.50 ~ -10.00D (in 0.50)		
Cyl			-0.75D	-1.25D	-1.75D
Axis			180° 170° 160° 90° 20° 10°	180° 170° 160° 90° 20° 10°	180° 170° 160° 20° 10° 0°



premiO

Recommended wearing time

Daily wear or flexible wear (up to 6 night extended wear) in 2 week replacement.

Lens Care

MeniCare™ Soft Multipurpose solution for cleaning, disinfection, storage and rinsing.
Innovative formulation: Comfortec™

- Efficient to remove lipids deposits.
- Less protein adhesion.
- High disinfection power with low cytotoxicity.

Club PremiO

All about Menicon PremiO lenses, e-memo service:
www.meniconpremio.com

Packaging

- 6 lens blisters per pack.
- Menicon heart shaped blister.
- 6 month PremiO and PremiO Toric Pack: 4 boxes x 6 lenses, + 3 x 360 ml MeniCare™ Soft, + 3 lens cases.



The 2 week silicone hydrogel lens

Made with Passion

