

# Menicon Soft MA


**Material : Govafilcon A**

Copolymere d'hydroxyethyl methacrylate  
 37% water content, non ionic  
 Blue handling tint  
 Dk:  $10 \times 10^{-11}$  ( Fatt )

**Sphero-aspheric design**

Soft hydrophilic contact lens

Dk/tc:  $11,1 \times 10$  (Fatt)  
 ( $e_c$  : 0.09 mm pour -3.00 D).



**Range**

Parameters	$r_0$ (mm)				Spherical power (D)	Diameter (mm)
		7.80	8.10	8.40		
step	0.30				0.25	-

**Fitting rules**

**In first intention  $r_0 = 8.40$  mm and  $\varnothing_T = 13.50$**   
 Diameter 13.00 mm for small cornea



84 M :  $r_0 = 8.40 - \varnothing_T = 13.50$   
 84 :  $r_0 = 8.40 - \varnothing_T = 13.00$

You can read the Menicon logo and BC number  
 Photochemical labeling.

**Indications**

- Soft hydrophilic lens, for myopia and hypermetropia correction.
- **Daily wear.**
- **Recommended lens replacement < 1 year.**

**Lens care**

**MeniCare™Soft**, Multipurpose solution for cleaning, disinfection, storage and rinsing.



**Packaging**

A lens vial containing sterile saline solution.