



Menicon PremiO

Material : Asmofilcon A

Silicone Hydrogel (Group V ISO18369-1).
 Water content: 40%.
 Dk : 129×10^{-11} (cm²/sec) mL O₂/(mL . mmHg)).
 Spectral and luminous transmittance: 98%.
 Refractive index: 27° (Captive bubble method).
 Blue handling tint.
 Nanogloss™ surface treatment.



Sphero-aspheric design

Dk/e : 161×10^{-9} (e_c : 0.08 mm for -3.00 D).
 Maximum corneal oxygenation across all parts of the lens and over the entire power range.
 Comfortable thin design with semi-round edges.

Range

Parameters	r ₀ (mm)	Spherical power (D)	Diameter (mm)
Standard	8.60 et 8.30	+ 6.00 to - 6.50 per 0.25 D - 6.50 to - 13.00 per 0.50 D	14.00

Fitting rules

BC : 8.60 mm as first choice.

Lens mark with the number 6 or 3 to read inside the lens.



Indications

- For myopia and hypermetropia.
- Daily wear or flexible wear (up to 6 night extended wear).
- 2 week replacement.

Lens care

MeniCare Soft™ Multipurpose solution for cleaning, disinfection, storage and rinsing.



Packaging

6 lens blisters per pack.
 Menicon heart shaped blister

