

Questions on age-related far-sightedness (presbyopia)

? You already wear GP lenses

Then you are already convinced of the advantages offered by this type of rigid O2-permeable lens (**comfort, safety, quality of vision etc.**).

Don't hesitate: contact your specialist to find out which lenses to try out and define their new characteristics, including correcting your far-sightedness.

? You have never worn contact lenses

If there is no age for starting to wear contact lenses, so why not try the new **Menicon Z multifocal lenses** ?

Consult your eye care specialist, who will explain all of the possibilities on offer and help you discover Menicon lenses. Wearing ultra-permeable contact lenses made out of Z material is choosing **the safest option** for the health of your eyes.

? Your eyes are becoming more sensitive

As you get older, it's normal to have fewer tears. Z ultra-permeable contact lenses do not dehydrate and let your eyes breathe: they are **first-line treatment** contact lenses that favour **long-term tolerance**. You can wear this type of contact lens even if you have few tears or feel as though you have dry eyes in air conditioned rooms or when in front of a screen. For "sensitized" cornea, it is the No. 1 and most reassuring option.

? You are already short-sighted or far-sighted and/or astigmatic

Thanks to the exceptional properties of Z material, Menicon offers a range of varied geometries that can adapt to different profiles and requirements : **from the simplest to the most complex**. Quality of vision is unequalled and more stable over time. Far-sightedness is an additional sight



defect that will be compensated at the same time as your ametropia with Z multifocal contact lenses.

? Handling the lenses

Z ultra-permeable contact lenses are **easier to handle**. Inserting and removing them involve **simple gestures** that subsequently become a mere reflex.

Wearing your lenses continuously for up to 30 days is one solution if you want to completely forget about your sight defect and not have to handle your contact lenses every day. This option is possible with Z material (Menifocal Z etc.) and may be envisaged with the consent of your eye care specialist. For the safety of your eyes, it is important to fully comply with all advice and instructions you are given.

Properly maintaining your lenses is important ...

- ✓ To keep your eyes in good health.
- ✓ To guarantee optimum quality of vision
- ✓ To optimize your comfort and the lifetime of your contact lenses.

It's simple with :



www.menicon.com



Your close-up vision is diminishing with age ?

You are becoming far-sighted !

Discover the new Z multifocal contact lenses
Discover the new Z multifocal contact lenses
Discover the new Z multifocal contact lenses
Discover the new Z multifocal contact lenses



Menicon contact lenses for clear distant and near vision



As soon as the first signs of age-related far-sightedness appear (blurred close up vision, particularly under subdued light, reading from further away etc.), consider Menicon Z multifocal contact lenses.

Since your far-sightedness changes over time, Menicon Z multifocal contact lenses change as well, guaranteeing you optimum quality of vision.

A range of Z multifocal contact lenses for all types of short-sightedness and for all stages of age-related far-sightedness :

- You have never worn either glasses or contact lenses
- You already wear glasses for far-sightedness
- You are already short-sighted or far-sighted and/or astigmatic.

A contact lens solution especially for you !

After an eye examination and listening to your near and distant vision requirements, your specialist will determine exactly what type of lens is most suitable for your eyes and lifestyle. He or she will define the characteristics of your contact lenses by determining the addition required for your far-sightedness (if you already wear glasses, it may be different).

Your specialist will tell you how to wear and maintain your lenses and how frequently you have to check them. You will very quickly learn how to insert and remove your contact lenses by wearing them progressively.

All Menicon contact lenses are made of **Z material, ultra-permeable to oxygen**, allowing them to be worn daily in all manners and/or continuously for up to 30 days in complete safety.



Each contact lens has a very specific near view and distant view zone. Several choices of geometry are proposed according to two vision principles: Alternate vision or simultaneous vision.

Alternate vision...

You decide whether you want to see close-up or far-off; it all depends on where you direct your gaze.

While looking straight ahead



Distance vision



While looking downwards (without lowering your head)




Near vision



Based on the principle of alternate vision, two types of contact lens are available: Menifocal Z and Menicon Z Executive.


Concentric contact lenses :

Menifocal Z



Segmented contact lenses :

Menicon Z Executive



Simultaneous vision...


All images reach your retina at the same time, whatever the distance and the object you are looking at: it is your brain that spontaneously sorts out these images, only conserving clear distant or near images.



Distance vision





Near vision



Menicon progressive contact lenses, based on simultaneous vision, are available in several designs even for astigmatism.

Progressive contact lenses :

Menicon Z Progressive / Z Progressive Near



- Distance vision
- Near vision
- Progression zone
- Transition zone

Forget your far-sightedness... Try Z multifocal contact lenses.